Light and Healthy

Light is a prerequisite for all living healthy growth

It's not only plants who need light . all creatures can not leave light . As early as the 2nd century we had the records of sunlight to treat disease .(solar therapy: sunlight, enhanced physique, prevent rickets).Artificial light therapy began at the end of the 18th century, to the middle of 19th century, visible light, infrared, ultraviolet, such as quickly formed, then the clinical treatment of widely in various fields of application and development.

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Introduction

LED optical mask made of high purity, high power density of red light, blue light to illuminate of skin, can change the cellular structure, killing bacteria, to provide a suitable environment for new cells, enhance new collagen's elastin, promote cell growth;Complex inflammatory acne skin aging, and alleviate the sun burns skin, and not do harming the skin, can whiten skin, improve skin elasticity.

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
Product principle - Wide spectrum of strong pulse photon technology

Specific wavelength pulse light energy can be through the epidermis, selective effects on the skin pigment, open blood vessels, hair follicles, atrophy, decomposition and absorption.At the same time also can make the thickening collagen tissue, generally after 15days to achieve whitening and anti-wrinkle treatment.Photon tender skin belongs to Non-ablative treatment and comfare tender skin to before (chemical peeling, skinned paste), it won't cause adverse effects and complications, safe and reliable.

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
Introduction

Adopting the photoradiation theory, LED Mask can make the user maintain the facial skin safely and easily with intimate design. It can promote the growth of skin collagen to achieve greater depth cosmetology of skin. It is the phototherapy you can enjoy any time as you want In addition, comparing with the traditional facial mask, LED Mask is more convenient, safer, and can be reused and has the advantage of environmental protection and others. LED Mask delivers voluminous LED red light energy to skin during the use, to process the crease from greater depth. It can not only flat the skin, but also prevent the growth of intradermal indigenous bacteria. You only need wear LED Mask 25 minutes every day, then LED red light will improve your skin gradually.

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



Product Benefits

Juvenility: Enhance the generation of collagen, reduce microgroove and wrinkles, and lighten skin color.

Compactness: Tighten the pores, improve the overall skin brightness, make the skin smooth and supple.

Rejuvenation: Fade freckles and opcs-redness, reduce acne and inflammation, accelerate the improvement of skin and the recovery after sunburn.

Moisturizing: Promote the circulation of face blood, lock up the skin aqua.

Vitality: Provide the cell activation energy, accelerate skin metabolism rate, thus to make the skin tender like a new-born baby.

Reduce acne: Restrain the inflamed phenomenon, accelerate the improvement of defective skin and repair of burn skin, and firm skin.

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LED蓝光面膜 / LEDブルーライトマスク / LED-Blaulichtmaske

Introduction

Blue light wavelength has the function for inhibiting melatonin.Light can affect the body's hormone secretion, heart rate, alertness, sleep quality, body and gene expression, blue-ray elevated body temperature, speed up the heart rate and reduce the role of sleepiness. Blue light can inhibit sleep related delta brainwave, alpha brain-wave increased significantly related to alertness.This means that you can use the short wave-length light blue as lethargy of control measures.LED Blue-ray can have antiseptic and anti-inflammatory effects to improve the surface of the skin, have very good effect to treat acne and rosacea.Blue-ray is specifically kill cause acne acne bacillus, just need exposure 15-25 minutes on the skin, only for acne, light is very gentle, sensitive skin can be done, without any side effects.For the red and blue light acne r side effects, there is no need to worry about it.

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



LED Blue Mask
LED蓝光面膜 / LEDブルーライトマスク / LED-Blaulichtmaske

Product Benefits

Bacteriostat: Repairing of acne damage and inhibition of wound bacteria.

Smoothing: Reduce the odds of scar caused by acne wound.

Anti anaphylaxis: Slow down the skin inflammation.

Glossiness: Improve the glossiness of skin.

Health: Reduce somnolence, increase vitality.

Conditioning: Adjust the biological clock (sleep disorders, time difference, Seasonal obstacle).

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Russia Certification

Qualifications & Patent

We have patent for LED mask , access to qualified EU market, advanced production systems, improved the quality of management

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Product quality inspection report

Our LED mask pass the test by China Quality Supervision and Inspection Center, Product safety non-toxic meet, safe household products, and skin and hair care appliances - Particular requirements, exceeding the national standard product requirements .

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
Method one: Connect to the computer
Take out the accessory cable of LED Driver, insert MICRO connector into DC-IN interface of LED Driver, and then connect the other end USB with the computer USB.
Method two: USB adaptor
Take out the accessory cable of LED Driver, insert MICRO connector into DC-IN interface of LED Driver, and then connect the other end USB with USB of DC 5V LED Driver.

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Matters Needing Attention

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Notices

- 1.Pregnant and minors are not suggested to use.
- 2.Part wavelength and frequency may cause epileptic disease, please do not use.
- 3.If suffering from thyropathy or taking thyropathy related medicines, please do not use.
- 4.If allergic to light or taking light-allergic medicines, for example,tetracyline, please do not use.
- 5.Try to avoid looking in the eye of the lighting area when use.

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Mask Maintenance

- 1.Use clear dry soft cotton cloth to wipe the surface gently in daily maintenance.
- 2.Do not use any solvent to clean the mask.
- 3.Do not disassemble the mask.
- 4.Please do not crush, throw, knock the mask and the LED Driver, keep LED Driver and lighting part dry.

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

Light and fashion amusement

Lights deductive fashion rich entertainment life

With the improvement of people's living standard, people's life quality requirements also more and more high, from the family decoration, household items to the leisure entertainment, people begin to pay attention to the elements of fashion, now the most popular in life is the science and technology elements, and gorgeous lighting is the most fashionable way to display.

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Introduction

With the development of science and technology in recent years, people's life is full of colour, even a restaurant or bar, the simplest menu became fashionable. LED Winelist is the latest development of fashion product in our company, its high brightness and brightness even, the built-in sustainable charging LED Driver, change picture simply, product appearance fashion contracted, unmatched by similar products on the market.

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
Introduction

Ice bucket, just as its name implies, it is used for cooling wine. When temperature is above the best temperature for drinking, ice bucket can cool it down in a few minutes. Company's latest LED ice bucket, suitable for bars and other entertainment place. Product appearance vogue, structure firm, switching LED Driver can emit colorful lights, beautiful and striking.

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
Function

Use RGB strip, brilliant color, built-in lithium battery for easy movement, many of output interfaces, can supply the power to the various equipments, built-in wireless transmitter, can charging for the mobile devices with wireless charging. Suitable for the use of KTV bar and other entertainment place.

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

Function

Use high light led strip, high brightness, built-in lithium battery for easy movement, many of output interfaces, can supply the power to the various equipments, the charger built-in the flexible plug, take out the plug from the box in the bottom of desk and connect to the LED Driver when charging and put it back to the bottom of the desk when not charging. suitable for the use of KTV bar and other entertainment place.

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

Introduction

LED light glass, suitable for all kinds of entertainment place under the dark condition. It can shine bright lights immediately as long as you pour drinks. Cooperating with glasses of wine, LED light glass will let you feeling in the future. Colorful lighting glasses with the entertainment lighting effects, will give you a deep impression. And for security reasons, whole cup materials using PVC plastic material, to ensure that a drunken customers will not be broken glass and scratch the body.

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
Introduction

LED tray, used in bars, dance halls and other places of entertainment, when the waiter carrying light tray send wine to the guest, extraordinarily appear in the dim light, the light emitted by the waiter's hand tray, produced by the refraction of light though the wine bottle, wine bottle also shine at the same time. It's an essential product in the bar and ballroom. Products use mobile power (or batteries), high brightness LED chip, high quality light guide, to ensure that the products are highlight and using waterproof processing to prevent the liquid on the tray damage light-emitting circuit board.

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
Introduction

LED dance floor is made of digital controlled color luminous plate and toughened glass. It can carry people for dancing, walking and other activities. Widely used in stage, places of entertainment, commercial square, Pedestrian Street, etc. It also can be used for exhibition to attract the eyeball. Products through dynamic, digital control color and brightness, can create the decorative effect of static and dynamic and be increase entertainment atmosphere quickly.

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

Introduction

Perfect even lighting
Great reduction of transportation costs
The perfect combination of Neon chalk board and light panelFast start up, no flicker, no noise
Using RGB source,PMMA materials, more than 90% light output, even light,no glare

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